



Be Right™



Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

Nr produktu: 2096250

PLN Cena: Kontakt

Podanie daty wysyłki

Vials, Tubes, Bottles & Racks

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Parametry techniczne

Długość:	150 mm
Objętość/Rozmiar opakowania:	10 /pk
Pojemność:	Capacity to Neck 50 mL
Średnica:	25 mm