



Be Right™



MEL Portable colorimeter/pH/microbiology lab for potable water

Nr produktu: 251242

NIEAKTUALNA POZYCJA
Ten produkt nie jest już dostępny.

Potable water kits make field testing easy

Hach potable water labs provide all equipment necessary for microbiological testing of Total Coliforms and E. coli using P/A methods, including Hach's portable incubator. The kits also include apparatus and reagents for testing 8 key water quality parameters, along with a DR 900 multiparameter handheld colorimeter and pH & TDS POCKET PRO testers. Microbiological media has to be ordered separately.

Parametry techniczne

Czułość: CFU/mL